

## Protein phosphatase assay

This protocol describes the standard strategy for measuring Ser/Thr protein phosphatase (PPase) activity in our laboratory using an artificial substrate (ex. Fzy-S50), a recombinant protein kinase (ex. Cdk) and [ $\gamma$ - $^{32}$ P]-ATP. This protocol can be modified/utilized to measure various PPase activity of your interest by changing substrate and kinase.

This protocol consists from 3 parts.

- 1, Purification of substrate protein
- 2, Phosphorylation of substrate
- 3, Phosphatase assay in *Xenopus* egg extract

## Reference

- 1, Mochida and Hunt, 2007, *Nature*, 449(7160), 336-340, PMID 17882219
- 2, Mackintosh C. & Moorhead G. *Protein phosphorylation: A practical approach* (ed. Hardie D.G.) pp 153-181 (Oxford University Press.1999)

## **Purification of substrate protein**

This protocol uses an artificial substrate, which is composed of tag protein (such as maltose binding protein) and a short peptide sequence derived from native phospho-protein.

### **Artificial substrate**

Fuse 25a.a.-length phospho-peptide sequence to MBP or GST expression vector. Phosphorylation site should be at the center of this peptide. I strongly recommend you to start with 5 or more different peptide sequences. Because some sequences are not expressed in bacteria or are not phosphorylated efficiently *in vitro*. And I do NOT recommend putting more than 1 phosphorylation site in each substrate, as different sites could be dephosphorylated by different PPases.

### **Materials**

BL21(DE3) strain + pLysS +pMALc2-FzyS50)  
Amylose resin (New England Biolab)

### **Protein Expression**

1 L of LB media supplemented with ampicillin & chloramphenicol  
Add 1 mM IPTG to your culture and incubate it for another 3 hrs @37 °C

### **Protein purification**

Collect bacterial cells & wash cell pellet once in PBS  
Suspend them in 25 ml of Extraction buffer +lysozyme  
Freeze this solution using liquid nitrogen  
Thaw cell pellet and incubate it in 37°C water bath for 15 minutes  
Cool it down on ice  
Sonicate it for 20-30 seconds 3-4 times (until no viscosity)  
Centrifuge @10 krpm for 15 minutes @4 °C  
Apply the supernatant to 5 ml of amylose resin  
Rotate the mixture for 0.5 hr @4°C  
Wash the beads 3 times with the Extraction buffer  
Add 5 ml of Elution buffer to the beads and rotate for 10 minutes @R.T.  
Retrieve eluate  
Repeat elution twice more  
Concentrate 15 ml of eluate (5 ml x3) into 1-2 ml using a centrifugal filter unit  
Dialyze it 3 times against Elution buffer  
Measure protein concentration  
Stock it in 1-mg aliquots @-20°C  
You would get 10-20 mg substrate from 1L bacterial culture

**Extraction buffer** : 20 mM Tris-Cl pH7.5, 150 mM NaCl, 1 mM DTT, 1 mM PMSF, 2 mM Benzamidine, 1 mM EDTA, 1 mM EGTA, 0.1% Tween-20

**Elution buffer** : 20 mM Tris-HCl pH7.5, 150 mM NaCl, 10 mM maltose

## Phosphorylation of substrate

### Materials

1 mg of substrate protein  
Active recombinant kinase  
[ $\gamma$ - $^{32}\text{P}$ ]-ATP

### Phosphorylation reaction

[Reaction mix]	
1 mg	Substrate
1-5 micro-g	Kinase (depends on how active your kinase is)
5 micro-L	[ $\gamma$ - $^{32}\text{P}$ ]-ATP
40 micro-L	5x kinase buffer
To 200 micro-L	Distilled/deionized water

Incubate your reaction mixture at 37°C overnight (longer than 12 hrs)  
Add 200  $\mu\text{l}$  of amylose resin (New England Biolab) to the reaction mix  
Rotate for 15 minutes  
Wash beads 5 times in Wash buffer  
Add 200 micro-L of Elution buffer to your beads  
Rotate for 10 minutes  
Retrieve 1<sup>st</sup> eluate to a new tube  
Repeat elution twice more, resulting in 600 micro-L of eluate  
Remove free ATP by using a filter unit and Wash buffer  
(until radio activities of flowthrough does not change)  
Concentrate substrate solution into 100 micro-L  
Measure its radioactivity using a liquid scintillation counter  
If successful, your preparation would have 100 ~200 kcpm/micro-L  
Freeze in 10-20 micro-L aliquots @-20°C

**5x CDK kinase buffer** : 0.1 M Hepes-Na pH7.8, 50 mM  $\text{MgCl}_2$ , 75 mM KCl, 5 mM EGTA, 25 mM NaF, 0.1 M  $\beta$ -glycero  $\text{PO}_4$ , 50 micro-M ATP, 5 mM DTT (Very low concentration of cold ATP is to get highly radio-active substrate)

[ $\gamma$ - $^{32}\text{P}$ ]-ATP : we use PerkinElmer NEG502A, 3000 Ci/mmol, 10 mCi/mL

**Wash buffer** : 20 mM Tris-Cl pH7.5, 150 mM NaCl, 0.05% Tween-20

**Elution buffer** : Wash buffer +10 mM maltose

**If your [ $\gamma$ - $^{32}\text{P}$ ]-ATP is old**, do increase its amount used accordingly.

## Phosphatase assay in *Xenopus* egg extract

### Materials

*Xenopus* egg extract (or your PPase source)  
Phosphorylated substrate

### Reaction pre-mix for 10 assays

200 kcpm            PPase substrate  
1 micro-L          Okadaic acid\* (100 micro-M in DMSO) or DMSO  
To 10 micro-L    Buffer (20 mM Tris-Cl pH7.5, 150 mM NaCl)

### PPase reaction

Mix egg extract (3 micro-L) and reaction pre-mix (1 micro-L)  
Incubate it @23°C for 4 minutes  
Add 20 micro-L of 10 % trichloroacetic acid to your reaction  
Vortex it for 5 seconds (you see white protein precipitate)  
Incubate it on ice for longer than 10 minutes  
(You can leave them on ice until all samples become ready)  
Centrifuge it for 5 minutes @13 krpm @4°C  
Retrieve all supernatant into a new tube, while trying not to take precipitate  
Centrifuge it for another 10 minutes @13 krpm @4°C  
Retrieve 15 micro-L of its supernatant into a new tube  
Add 20 micro-L of 5% (w/v) ammonium molybdate\*\* (in 0.5M H<sub>2</sub>SO<sub>4</sub>)  
Add 80 micro-L of organic solvent mix\*\*\*  
(2-methyl-1-propanol : heptane = 50 : 50, water-saturated)  
Vortex it vigorously for 30 seconds  
Centrifuge it for 1 minute @13 krpm  
Take 50 micro-L of organic (upper) phase to a scintillation tube with appropriate scintillant (we use Ecoscint A)  
Shake it well by hand  
Measure radioactivity using a scintillation counter

\* 2.5 micro-M of okadaic acid can suppress most of PP2A activity in *Xenopus* egg extract. You need very high concentration of this kind of inhibitor, because, in crude cell extracts, PPases are expressed at concentrations much higher than their reported IC<sub>50</sub> values.

\*\* I use SIGMA #A7302, but products from other suppliers would work too.

\*\*\* Molybdate/organic solvent extraction method is introduced in detail in Mackintosh C. & Moorhead G. *Protein phosphorylation: A practical approach* (ed. Hardie D. G.) pp 153-181 (Oxford University Press.1999).

**IMPORTANT !!!** It is crucial to keep protein concentration as high as *in vivo*. Because simple dilution of cell extract disrupt physiological regulation of PPases by unknown mechanism.